

Title:	Cooler Compressor Conditions			ID:
				127
19.03.13	Response:	Model:	Author:	
		-	CMA	

STC Cooler Unit Compressor Activate Conditions

Firmware

Version V3.0.00

Main Unit (3000=1)

1 1st Stage R404A

Contact	Type	Status	Comment
5003	Relay	1	<i>Delay T41 after power on (10s)</i>
1403	Flip-Flop	1	<i>Triggered by 5003 if not in service mode</i>
013	Relay	0	<i>Manual Switch</i>
006	Relay	1	<i>1st stage Compressor Over-Temperature Switch</i>
007	Relay	0	<i>1st stage Compressor Low-Pressure Switch</i>
3408	Flip-Flop	0	<i>Triggered by 008 1st stage Compressor Over-Pressure Switch, Reset by 1900</i>
014	Relay	1	<i>Power Voltage ok</i>
503	Relay	1	<i>Output to 1st stage Compressor Power Relay</i>

2 2nd Stage R508B

Contact	Type	Status	Comment
6500	Relay	1	<i>Stepper Valve successfully initialized</i>
T4	Relay	1	<i>Timer T004 delay 1st stage running, triggered by 004, 300s</i>
1703	Relay	1	<i>Temperature refrigerant < -10°C, DM984<DM894, delayed by T22</i>
T45	Relay	1	<i>Delay 10s</i>
1709	Flip-Flop	0	<i>2nd Stage Over-Temperature, triggered by DM983>#1100 and delayed by T29 (3600s), reset by power-cycle</i>
5004	Relay	1	
1404	Flip-Flop	1	<i>Triggered by 5004 if not in service mode</i>
013	Relay	0	<i>Manual Switch</i>
010	Relay	1	<i>2nd stage Compressor Over-Temperature Switch</i>
011	Relay	0	<i>2nd stage Compressor Low-Pressure Switch</i>
012	Relay	0	<i>2nd stage Compressor Over-Pressure Switch</i>
014	Relay	1	<i>Power Voltage ok</i>
504	Relay	1	<i>Output to 2nd stage Compressor Power Relay</i>

Auxiliary Unit (3001=1)

1 1st Stage R404A

Contact	Type	Status	Comment
5003	Relay	1	<i>Manual switch 000=1 or Error input 003=0 or Defrost Feed-Back 009=1 delayed by T41 (10s)</i>
1403	Flip-Flop	1	<i>Triggered by 5003 if not in service mode</i>
013	Relay	0	<i>Manual Switch</i>
006	Relay	1	<i>1st stage Compressor Over-Temperature Switch</i>
007	Relay	0	<i>1st stage Compressor Low-Pressure Switch</i>
3408	Flip-Flop	0	<i>Triggered by 008 1st stage Compressor Over-Pressure Switch, Reset by 1900</i>
014	Relay	1	<i>Power Voltage ok</i>
503	Relay	1	<i>Output to 1st stage Compressor Power Relay</i>

2 2nd Stage R508B

Contact	Type	Status	Comment
6500	Relay	1	<i>Stepper Valve successfully initialized</i>
T4	Relay	1	<i>Timer T004 delay 1st stage running, triggered by 004, 300s</i>
1703	Relay	1	<i>Temperature refrigerant < -10°C, DM984<DM894, delayed by T22</i>
T45	Relay	1	<i>Delay 10s</i>
1709	Flip-Flop	0	<i>2nd Stage Over-Temperature, triggered by DM983>#1100 and delayed by T29 (3600s)</i>
5004	Relay	1	
1404	Flip-Flop	1	<i>Triggered by 5004 if not in service mode</i>
013	Relay	0	<i>Manual Switch</i>
010	Relay	1	<i>2nd stage Compressor Over-Temperature Switch</i>
011	Relay	0	<i>2nd stage Compressor Low-Pressure Switch</i>
012	Relay	0	<i>2nd stage Compressor Over-Pressure Switch</i>
014	Relay	1	<i>Power Voltage ok</i>
504	Relay	1	<i>Output to 2nd stage Compressor Power Relay</i>