

Title:	Glass Door lose	ID:	
		0116	
Date in:	Response:	Model:	Author:
2003-06-13	2003-06-13	-	CMa

Q:

I'm not sure what type incubator they have. The FSE was at this site as well today for another issue. When there the glass door of the incubator feel off. Apparently the epoxy that mounts the hinges to the glass has come loose. This seems like a much better solution to mount the door. Would it be possible to get them the new door? If not, what can we do about the existing door and getting it fixed?

A:

This must be a very old system. We changed the bonding to screws. It would be helpful if you could send us the dimensions of the glass door. We can send you a replacement door.