

Title:	Different Plate Size per Stacker or Carrousel	ID:	
		<i>0081</i>	
Date in:	Response:	Model:	Author:
<i>2003-01-21</i>	<i>2003-01-21</i>	-	<i>CMa</i>

Q:

Dave's customer is asking if they can mix and match plate sizes within an incubator.

Is it possible to have both standard and deepwell plates in the same incubator, or at least some racks of each?

A:

Yes, you can.

You must set the pitch and the max. number of plates per stacker correctly. Some people send the pitch (WR DM23 pp) on every access.