

Title:	Liconel Appearance	ID:	
		0334	
Date in:	Response:	Model:	Author:
2008-05-23	2008-05-23	-	CMa

Q:

As far as we watch the chamber, it looks as standard stainless chamber, not Liconel. Does Liconel look like standard chamber ?
How he can judge whether it is standard or Liconel ?
He needs to confirm as soon as possible. Please get an information from Liconic AG, whether this unit was built with Liconel.

A:

Yes Liconel looks like stainless steel but it is slightly darker and less shiny. If you send me the serial number I can verify.