

Title:	Shaking of Yeast	ID:	
		0345	
Date in:	Response:	Model:	Author:
2009-11-23	2009-11-23		CMa

Q:


Can a 1mm Shaker keep yeast floating?

A:

The Pictures below show a Test plate with yeast. First picture shows the plate before treatment. The second picture shows the plate after less than two minutes in the 1mm Shaker at 1200 rpm.



Picture 1: Plate before treatment in the Shaker



Picture 2: Plate after less than 2 minutes in the 1mm Shaker at 1200 rpm