

Title:	Why use Bar Code Reader	ID:	
		0355	
Date in:	Response:	Model:	Author:
2010-04-21	2010-04-21		CMa

Q:

Why Bar Code Reader?

A:

a) BCR is typically used for

- Plate inventory
- BC counter-check

having two or one BCR will depend on where plates are typically fed into the system. If the user plan to manual put plates into the incubator and in the plate hotel, we would advise to have two BCR's. We typically would monitor the user doors. If we register a manual door access we assume that the system's plate content may have changed and we perform a complete inventory of the system's plate content.

b) This is a stand-alone application that has no external housing around it. We recommend a complete housing (not only half housing - like our old MPH that is now distributed by Thermo). In addition, the full housing provides better user protection and allows monitoring the user door and therefore detecting manual user accesses.