

<b>Title:</b>	<b>Maximum Power Consumption STX44-HRxx</b>			<b>ID:</b>
				0380
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2004-04-27	2004-04-27	STX44	C <i>Ma</i>	

**Q:**

What is the maximum power consumption of a STX44 HRBT if we start from ambient temperature and we ask 4°C (91180900). I could not find anywhere in the manual.

**A:**

~600W. Note, there is a high peak on start-up of compressor. We use D-type circuit breakers.