

Title:	Temperature-RH Sensor Probe			ID:
				133
07.12.13	Response:	Model:	Author:	
		-	CMA	

Temperature / RH Sensor

There are two types of Temperature / RH Sensors (Probes) built in to climate controlled units. The old type sensor has been discontinued. The transition took place in 2011.



The probes are not interchangeable and have

- ! Different Connector
- ! Different Sensor Cable
- ! Different Display (DDC)

In case of any sensor failure the sensor must be replaced by the same type sensor. A limited number of old type sensors are on stock at Liconic.

The calibration of the new type sensor requires a special communication cable.