

Accessories – Transfer Stations

BT-Attachment Kit



This transfer nesting station is our standard design and works on all of the LiCONiC Products. The kit can be used for customers who have changed the design of their work-station and need a standard solution for their next design.

The components used are the most reliable and highest quality products available. All of our manufactured components are electro-coated for the prevention of corrosion and contamination from external conditions.

Specifications

Available for	All Bench Top and Stand Alone units of the STX Series except STX44-xxSA All remaining automated incubator series except STT1000-Kiwi. All plate hotel series
----------------------	--

Design	Fixture for Transferstations
---------------	------------------------------