

## Accessories

### Cassettes

All members of the StoreX series store the plates in cassettes that can easily be handled through the bis user door by the operator.

Cassettes are available for various plate dimensions.

Even customized sizes are available at only slightly extra cost. All-metal cassettes are also available.

#### Function

For the storage of the microtiter plates a StoreX 44 system is equipped with 2 and a StoreX 220 system is equipped with 9 stainless steel magazines. The number of levels within the magazine is depending on the type of micro plate.

#### Main features

- stable construction
- stainless steel
- can be taken out of the system for loading, unloading or cleaning
- can be exchanged between different systems
- storage of microtiterplates with lids possible

#### Specifications

Plate Type	Plate Height	Pitch	# Plates Cassette		
			505 mm	645 mm *	1210 *
NTP (Nanotiter Plates)	5 mm	11	42	55	111
	11 mm	17	28	37	72
1536	12 mm	18	27	35	68
MTP (Microtiter Plates)	17 mm	23	22	28	53
	22 mm	28	17	23	43
	23 mm	29	17	22	42
	24 mm	30	17	21	41
	27 mm	33	15	19	37
DWP (Deep Well Plates)	44 mm	50	10	13	25
	53 mm	59	8	10	21
DiTi200	66 mm	72	7	8	17
DiTi1000	104 mm	110	4	5	11

\* can only be used with STX500 or STX1000 bi-carousels