

## Accessories – Transfer Stations

### Dual Transfer Station



The dual slide-station was designed for use by non-articulating robotic arms which need additional clearance from the storage units. With its precise stopping point it ensures that the devices need to be taught only once for accurate load and retrieval of the storage plates. The dual slide station can be used to place an STR Series (Periscope) product between two existing instruments. Plates from an STR Unit can be fed to both the external instruments simultaneously.

The components used are the most reliable and of the highest quality available in the market. All manufactured components are electro-coated to protect against corrosion and contamination from external environmental conditions.

#### Order information

STR44, STR240	<u>Order Nr.</u> 9143 01 58
---------------	--------------------------------

#### Specifications

<b>Plate Positions</b>	2
<b>Sensor</b>	Plate presence sensor
<b>Layout</b>	Straight on
<b>Displacement</b>	Perpendicular to and from the Unit
<b>Application</b>	Installation of a STR between two instruments