

Accessories – Transfer Stations

Flask Transfer Station

The flask transfer station allows the transport of “robo-flasks” and 175 ml T Flasks to your tissue culture work-station. It is the only design for this type of transport and work on the market. Automating your tissue culture process is now possible. The components used are the most reliable and highest quality products available. All of our manufactured components are electro-coated for the prevention of corrosion and contamination from external conditions. The belt driven system ensures smooth transfer and accuracy of delivery.

Specifications

Plate Positions	1
Sensor	Plate presence sensor
Layout	Perpendicular rail
Movement	Perpendicular to / from Unit
Special	Only for Flask Incubator