

Accessories – Transfer Stations

Removable Transfer Station



The Removable Transfer Station is available only for LiCONiC's STR Series. With its precise mounting point it ensures that the devices need to be taught only once for accurate load and retrieval of the storage plates. This Transferstation can be detached from the storage unit for easy service and maintenance access especially for under the bench units.

The components used are the most reliable and of the highest quality available in the market. All manufactured components are electro-coated to protect against corrosion and contamination from external environmental conditions.

Order information

Removable Transfer Station	<u>Order Nr.</u> 9143 01 72
----------------------------	--------------------------------

Specifications

Plate Positions	1
Sensor	Plate presence sensor
Layout	Nest
Displacement	None
Application	Easy maintenance access for STR units