

Accessories

StoreX Panel

The 'StoreX Panel' is a software package designed to control all types of StoreX incubators.

It allows to perform:

- Climate control, including climate logging and drawing of Graphs
- Plate handling operations
- Incubator gate and carousel control
- StoreX production tests
- Executing low-level direct commands

Specifications

Functions:	Climate control & logging Plate Handling operations Unit carousel and gate control Test routines Low level commands terminal
Climate:	Logging: Data stored as TXT Files Graphs: Data display Climate-Window: Read and Set Climate
Plate handling:	Drag & Drop or popup menu controlled Buttons for Gate and Carrousel control
Test routines:	Perform sequential Test: - Quick: first level – last level - All Slots Perform random access Test Output stored as TXT File
Low level commands:	Terminal window Commands history