

## Accessories

### Manual Transfer Station Nest

Special nest for manual single plate feeding operation. The Manual Transfer Station Nest has extra wide alignment cones for safe and simple placement of a plate into the transfer station nest. The Manual Transfer Station Nest includes a plate detection sensor for sensing plates placed in the Transfer Station Nest.



#### Specifications

<b>Plate Positions</b>	1
<b>Sensor</b>	Plate detection Sensor
<b>Layout</b>	Nest
<b>Movement</b>	None
<b>Plate format</b>	SBS
<b>Catching range</b>	+/- 6mm
<b>Material</b>	Aluminium, POM