

Accessories

Transfer Nest

Universal Transfernest for one plate and common interface between Liconic Instrument and surrounding system. Transfernest will accept SBS format plates according.

- ANSI/SBS 1-2004 Footprint Dimensions for Microplates
- ANSI/SBS 3-2004 Bottom Outside Flange Dimensions for Microplates

Specially developed alignment comes with optimised sliding properties will align plates for accurate transfer positioning. The Transfernest will correct for mis-alignment of over to +/- 2 mm.



Specifications

Plate Positions	1
Sensor	None
Layout	Nest
Movement	None
Special	