

Option

High Temperature

Description

The StoreX Series can be upgraded to operate at temperatures up to 70°C. The High Temperature option contains the markets most heat resistant material and at the same time withstand almost all organic and inorganic substances. This ensures long-term stability and will help your discoveries to be true and accurate.

For other higher temperature versions please contact our local representation.

Highlights of high temperature option

- Thermo-resistant material
- Chemically stable to probes
- High long term stability

Order information

<i>High temperature Option for ...</i>	<i>Order Nr.</i>
STX44, STR44	9118 11 02
STX110, STX220, STR240	9122 06 02
STX500, STX1000, STR602	9122 06 06
STF130	9141 02 01