

# Option

## StoreX Liconel

The New Liconel based StoreX Series is the markets most contamination resistant incubator family. Liconel technology combines de-contamination properties like copper with the advantage for extremely low aging effects like stainless steel. This insures long-term stability and will help your discoveries to be true and accurate.

### What is Liconel?

Liconel is an innovative material specially developed for the incubation chambers.

Liconel is a material that prevents unwanted microbial growth while maintaining protection of probes.

Liconel was specially developed for incubators with integrated handling.

### Highlights of Liconel

- Self – decontaminating
- Non – toxic to probes and user
- High long term stability
- Simple maintenance

### Order information

<i>StoreX Liconel for ...</i>	<u>Order Nr.</u>
STX40	9118 12 43
STX110	9131 05 11
STX220	9122 05 25
STX500	9132 06 20
STX1000	9136 05 11
STR44	9118 12 44
STR240	9134 01 97
STR602	9141 01 84

### Specifications

- Material: Alloy
- Principle of operation: Cu++ Ions
- Effect: Contamination free surfaces
- Application: Entire unit, Complete Chamber or only Critical parts