

Option

Plate Presence Sensor

General Information

The plate presence sensor has two functions:

The PPS allows the user to rapidly and conveniently take an inventory of the incubator or plate hotel and identify where plates are within the cassette format. This allows a quick and easy way to determine the available capacity of the unit to prevent any over-rides.

The plate presence sensor can also be used in conjunction with the barcode reader by identifying where the plates are within the units. This allows the barcode reader to create a much faster inventory by reading only those positions within the unit where plates are present saving valuable time.

Order information

	<u>Order Nr.</u>
PP-Sensor	9121 00 62

Specifications

- Use: Detection of Plates in Stacks
- Principle of operation: Reflective light barrier
- Application: Detection of occupied position
- Remarks: Standard component for Barcode reader option