

LPT220-EVO

Description


The LPT220-EVO is specially designed to work with Tecan platforms. It is the successor of the legacy LDR 189 plate capacity plate hotel.



Capacity

Plate Type	Plate Height	Capacity
NTP (Nanotiter Plates)	5 mm	420
	11 mm	280
1536	12 mm	270
MTP (Microtiter Plates)	17 mm	220
	22 mm	170
	23 mm	170
	24 mm	170
	27 mm	150
DWP (Deep Well Plates)	44 mm	100
	53 mm	80
DiTi200	66 mm	70
DiTi1000	104 mm	40

Dimensions



Width (W)	725 mm
Depth (D)	583 mm
Height (H)	876.5 mm
Transfer Height (T)	138.5 mm