

Option

LPR Housing



- + Enclosure for LPR Plate Hotels
- + large Front Door for user access
- + protection from unwanted access and light protection

This Option maybe retro-fit in the field

The LPR housing is designed to protect the complete robotics and prevent any operational fault due to external interference. This cover has a user door for an easy access to the stacks. The design of the cover is also suitable for light sensitive samples.

Order information

LPR240 Housing

Order Nr.
9143 01 36

Specifications

- Design: Full housing with door
- Protection: All robotic parts
- Plate Access: Plate transfer station and door for access to plate stacks