

## Option

### **LPX Cover Protection**



- + Enclosure for LPX Plate Hotels
- + large Front Door for user access
- + protection from unwanted access and light protection

This Option maybe retro-fit in the field.

The LPX cover protection is designed to protect the complete robotics and prevent any operational fault due to external interference. This cover has a user door for an easy access to the stacks. The design of the cover is also suitable for light sensitive samples.

#### **Order information**

	<u>Order Nr.</u>
LPX44 Cover Protection	9147 00 19
LPX220 Cover Protection	9144 00 56
LPX440 Cover Protection	9144 00 58

#### **Specifications**

- Design: Full housing with door
- Protection: All robotic parts
- Plate Access: Plate transfer station and door for access to plate stacks