

Option

Orbital Shaker

Agile, Popular & Innovative

For LPX44

General

Orbital shaker with integrated auto-positioning for handler access. All functions required for access implemented autonomously. Driver compatible with non-shaker unity, no special commands required for plate handling. All sub functions autonomously implemented.

Order information

<i>Orbital Shaking Option for LPX44</i>	<u>Order Nr.</u>
1mm Amplitude, 600 rpm max.	9147 00 16
1mm Amplitude, 1200 rpm max.	9147 00 17
3mm Amplitude, 600 rpm max.	9147 00 35
6mm Amplitude, 400 rpm max.	9147 00 36

Specifications

Type	Orbital	
Amplitude	1 mm	0-1200 rpm
	3 mm	0-600 rpm
	6 mm	0-400 rpm
Frequency	Software Adjustable	
Positioning	Integrated	
Handling Access Control	Automatic	
Manual Access	User Door	
Activate / Deactivate	Remote Control	

Option

Linear Shaker

Agile, Popular & Innovative

For LPX220, LPX440, LPR240

General

By accelerating and decelerating the carousel movement, a linear shaking function can be integrated in our carousel based units.

Via software set points the amplitude and the speed can be adjusted.

Order information

	<u>Order Nr.</u>
Linear Shaking Option	9122 05 15

Specifications

Type	Linear / Agitation
Amplitude	1-100 mm 12-300 rpm
Frequency	Software Adjustable
Positioning	Integrated
Handling Access	Automatic
Manual Access	User Door
Activate / Deactivate	Remote Control