

Option

Ultra High Speed

Function

The typical access time is 15 - 30 s. For HTS applications a high speed with shorter access time is offered.

Order information

<i>UHS for ...</i>	<u>Order Nr.</u>
LPX44	9111 10 00
LPX220	9111 10 00
LPX440	9111 06 04
LPR240	9111 10 00

Specifications

Accelerated handling operations

Average access time:

LPX44: 15 sec	LPX44 UHS: 12 sec
LPX220: 18 sec	LPX220 UHS: 12 sec
LPX440: 18 sec	LPX440 UHS: 14 sec
LPR44: 20 sec	LPR44 UHS: 18 sec
LPR240: 24 sec	LPR240 UHS: 19 sec