

# WDU

## Description

LiCONiC has created the perfect Waste Disposal Unit to protect workers and assays from contamination. The WDU can store up to 900 plates and lids, which can be retrieved from the unit for easy and safe disposal. The quick release cassette can be autoclaved or simply use the quick release to deposit plates gently into bio-contamination containers.




The usual method is to drop plates into a box which causes splashing and poses risk to end-users and airborne contaminants to assays. When safety is a concern, you should take all measures to ensure the safety and welfare of your employees and ensure accuracy of your assays.



## Capacity

| Plate Type         | Capacity |
|--------------------|----------|
| 384 – well plates  | 500      |
| 1536 – well plates | 750      |
| Lids               | 1200     |

## Dimensions



|                     |           |
|---------------------|-----------|
| Width (W)           | 715 mm    |
| Depth (D)           | 772 mm    |
| Height (H)          | 1325 * mm |
| Transfer Height (T) | 1030 mm   |

\* add 30-70mm for feet

## Specifications

### Weights and Dimensions

| Unit | Incl. Cassettes | Net Weight | Tare  | Gross weight | Package | Dimensions cm (LxWxH) |
|------|-----------------|------------|-------|--------------|---------|-----------------------|
| WDU  | None            | 135 kg     | 73 kg | 208 kg       | Crate   | 94 x 89 x 163         |